


SP-B Antibody / Surfactant protein B [clone RM370] (R20389)

Catalog No.	Formulation	Size
R20389-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM370
Purity	Protein A purified from animal origin-free supernatant
UniProt	P07988
Applications	Immunohistochemistry (FFPE) : 1:500-1:1000 Western Blot : 1:500-1:1000
Limitations	This recombinant SP-B antibody is available for research use only.



Western blot testing of human lung tissue lysate with recombinant SP-B antibody at 1:1000.



IHC staining of FFPE human lung cancer tissue with recombinant SP-B antibody at 1:1000.

Description

This antibody reacts to human SP-B (Surfactant protein B) and proSP-B.

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant SP-B antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A peptide corresponding to residues near C-terminus of human SP-B (Surfactant protein B) was used as the immunogen for the recombinant SP-B antibody.

Storage

Store the recombinant SP-B antibody at -20oC.