

SR-BI Antibody (R36278)

Catalog No.	Formulation	Size
R36278-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	949
Applications	Western Blot : 0.1-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This SR-BI antibody is available for research use only.



Western blot testing of human adrenal gland lysate with SR-BI antibody at 0.01ug/ml.
Predicted molecular weight: ~57/82kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: Scavenger receptor class B, member 1 (SCARB1); scavenger receptor B1

Application Notes

Optimal dilution of the SR-BI antibody should be determined by the researcher.

1. This SR-BI antibody is specific for isoform 1.

Immunogen

Amino acids TSAPKGSVLQEAK were used as the immunogen for this SR-BI antibody.

Storage

Aliquot and store the SR-BI antibody at -20oC.

References (1)