

SS18-SSX Antibody / Protein SSXT [clone SS18.SSX/12628R] (V5838)

Catalog No.	Formulation	Size
V5838-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	100 ug
V5838-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	20 ug
V5838SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG, kappa
Clone Name	SS18.SSX/12628R
Purity	Protein A affinity
UniProt	Q15532
Localization	Nucleus
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml Western Blot : 2-4ug/ml
Limitations	This SS18-SSX antibody is available for research use only.



Description

S18 is a protein that has been shown to be a part of the SWI/SNF complex. The SS18-SSX fusion proteins are a result of in-frame fusions that fuse the SS18 gene on chromosome 18 with X chromosome genes SSX1, SSX2, and to a lesser extent SSX4. Human synovial sarcoma (SS) accounts for 8-10% of all soft tissue malignancies and 95% of these

malignancies express the recurrent translocation of the SS18 gene on chromosome 18. The N-terminal SNH domain (SYT N-terminal homology domain) of the SS18 protein interacts with SWI/SNF chromatin remodeling complexes via the N terminal region of BRM and BRG1 subunits. Studies of the SS18-SSX fusion in SS suggest that endogenous SS18 competes with the mutant SS18-SSX fusion for occupancy in the SWI/SNF complexes resulting in the displacement of SNF5 (BAF47) subunit. Displacement of the SNF5 subunit results in altered function of the SWI/SNF complex that leads to deregulated expression of genes such as Sox2 in synovial sarcoma.

Application Notes

Optimal dilution of the SS18-SSX antibody should be determined by the researcher.

Immunogen

A recombinant fragment corresponding to the fusion site of human SS18-SSX protein was used as the immunogen for the SS18-SSX antibody.

Storage

Aliquot the SS18-SSX antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.