


STAT3 Antibody (F47582)

Catalog No.	Formulation	Size
F47582-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F47582-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Bovine, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P40763
Applications	Western Blot : 1:1000
Limitations	This STAT3 antibody is available for research use only.



STAT3 antibody western blot analysis in A375 lysate. Predicted/observed molecular weight ~88kDa.

Description

The protein encoded by this gene is a member of the STAT protein family. In response to cytokines and growth factors, STAT family members are phosphorylated by the receptor associated kinases, and then form homo- or heterodimers that translocate to the cell nucleus where they act as transcription activators. This protein is activated through phosphorylation in response to various cytokines and growth factors including IFNs, EGF, IL5, IL6, HGF, LIF and BMP2. This protein mediates the expression of a variety of genes in response to cell stimuli, and thus plays a key role in many cellular processes such as cell growth and apoptosis. The small GTPase Rac1 has been shown to bind and regulate the activity of this protein. PIAS3 protein is a specific inhibitor of this protein. Three alternatively spliced transcript variants encoding

distinct isoforms have been described.

Application Notes

Titration of the STAT3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 673-699 from the human protein was used as the immunogen for this STAT3 antibody.

Storage

Aliquot the STAT3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.