


STAT5A Antibody (R34705)

Catalog No.	Formulation	Size
R34705-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6776
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This STAT5A antibody is available for research use only.



Western blot testing of K562 lysate with STAT5A antibody at 0.1ug/ml. Predicted molecular weight: ~91kDa.

Description

Additional name(s) for this target protein: Signal transducer and activator of transcription 5A

Application Notes

Optimal dilution of the STAT5A antibody should be determined by the researcher.

Immunogen

Amino acids KYYTPVLAKAVD were used as the immunogen for this STAT5A antibody.

Storage

Aliquot and store the STAT5A antibody at -20oC.