


STMN1 Antibody [Stathmin 1] (R34545)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R34545-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|----------------------|---|
| Availability | 1-3 business days |
| Species Reactivity | Human, Mouse |
| Predicted Reactivity | Rat |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 3925 |
| Applications | Western Blot : 1-3ug/ml Enzyme ImmunoAssay : binds recombinant protein ELISA (peptide) LOD : 1:1000 |
| Limitations | This STMN1 antibody is available for research use only. |



Western blot testing of human testis lysate with STMN1 antibody at 1ug/ml. Predicted molecular weight: ~17kDa.

Description

Additional name(s) for this target protein: Stathmin 1

Application Notes

Optimal dilution of the STMN1 antibody should be determined by the researcher.

Immunogen

Amino acids KNKESKDPADETEAD were used as the immunogen for this STMN1 antibody.

Storage

Aliquot and store the STMN1 antibody at -20oC.