

Stomatin Antibody (R33847)


Catalog No.	Formulation	Size
R33847-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2040
Applications	Western Blot : 0.1-0.3ug/ml Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Stomatin antibody is available for research use only.



Western blot testing of human lung lysate with Stomatin antibody at 0.1ug/ml. Predicted molecular weight: ~32kDa.



IHC staining of FFPE human liver with Stomatin antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: STOM; BND7; EPB7; EPB72; erythrocyte band 7 integral membrane protein

Application Notes

Optimal dilution of the Stomatin antibody should be determined by the researcher.

Immunogen

Amino acids EKRHTRDSEAQR were used as the immunogen for this Stomatin antibody.

Storage

Aliquot and store the Stomatin antibody at -20oC.