


STRADB Antibody / Als2Cr2 / ILPIP (R33210)


Catalog No.	Formulation	Size
R33210-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	55437
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This STRADB antibody is available for research use only.



HEK293 overexpressing ILPIP/STRADB (lane 2) tested with STRADB antibody (lane 1-mock transfection). Predicted molecular weight: ~47/31kDa (isoforms alpha/beta).



IHC staining of FFPE human heart with STRAD8 antibody. HIER: steamed with pH6 citrate buffer, AP-staining.



Western blot testing of human heart lysate with STRADB antibody at 1ug/ml. Predicted molecular weight: ~47/31kDa (isoforms alpha/beta).

Description

Additional name(s) for this target protein: STE20-related kinase adaptor beta; ILP-interacting protein ILPIPA

Application Notes

Optimal dilution of the STRADB antibody should be determined by the researcher.

1. This STRADB antibody will detect both human isoforms of STRADB: isoform 1/alpha (418 AA) and 2/beta (280 AA).

Immunogen

Amino acids CDFPDEKDSYWEF were used as the immunogen for this STRADB antibody. (1)

Storage

Aliquot and store the STRADB antibody at -20oC.