


Survivin Antibody (R33364)

Catalog No.	Formulation	Size
R33364-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cat, Cow, Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	332
Applications	Western Blot : 0.01-0.5ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Survivin antibody is available for research use only.



Western blot testing of MOLT4 lysate with Survivin antibody at 0.01ug/ml. Predicted molecular weight ~16kDa.

Description

Additional name(s) for this target protein: Baculoviral IAP repeat-containing 5; BIRC5

Application Notes

Optimal dilution of the Survivin antibody should be determined by the researcher.

1. This Survivin antibody is specific to isoforms 1 & 3.

Immunogen

Amino acids CVRRAIEQLAAMD were used as the immunogen for this Survivin antibody.

Storage

Aliquot and store the Survivin antibody at -20oC.