


Synaptogyrin 4 Antibody (R34205)

Catalog No.	Formulation	Size
R34205-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23546
Applications	Western Blot : 0.01-0.1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Synaptogyrin 4 antibody is available for research use only.



Western blot testing of human brain lysate with Synaptogyrin 4 antibody at 0.01ug/ml.
Predicted molecular weight: ~26kDa.

Description

Additional name(s) for this target protein: SYNGR4

Application Notes

Optimal dilution of the Synaptogyrin 4 antibody should be determined by the researcher.

Immunogen

Amino acids HIPKSLQELANS were used as the immunogen for this Synaptogyrin 4 antibody.

Storage

Aliquot and store the Synaptogyrin 4 antibody at -20oC.