


Syntaxin 6 Antibody (R35196)


Catalog No.	Formulation	Size
R35196-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10228
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Syntaxin 6 antibody is available for research use only.



Western blot testing of human lymph node lysate with Syntaxin 6 antibody at 0.3ug/ml.
Predicted molecular weight: ~29kDa.



Western blot of HEK293 lysate overexpressing human STX6 probed with Syntaxin 6 antibody (mock transfection in lane 1). Predicted molecular weight: ~29kDa.

Description

Additional name(s) for this target protein: STX6

Application Notes

Optimal dilution of the Syntaxin 6 antibody should be determined by the researcher.

Immunogen

Amino acids QALAERKNRQA were used as the immunogen for this Syntaxin 6 antibody.

Storage

Aliquot and store the Syntaxin 6 antibody at -20oC.