


TAL1 Antibody / T-cell acute lymphocytic leukemia protein 1 [clone BTL73] (V8109)

Catalog No.	Formulation	Size
V8109-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	100 ug
V8109-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	20 ug
V8109SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1, kappa
Clone Name	BTL73
Purity	Protein G affinity chromatography
UniProt	P17542
Applications	Western Blot : 1-2ug/ml
Limitations	This TAL1 antibody is available for research use only.



SDS-PAGE analysis of purified, BSA-free TAL1 antibody (clone BTL73) as confirmation of integrity and purity.

Description

Activation of TAL1 characterizes up to 60% of cases of human T cell acute lymphoblastic leukemia, making it the most frequent gain-of-function mutation observed in this disorder. TAL1 (also designated SCL) is a serine phosphoprotein and

basic helix-loop-helix transcription factor known to regulate embryonic hematopoiesis. This transcription factor binds as a heterodimer with E2A and HEB/HTF4 to a nucleotide sequence motif termed the E-box. In addition, leukemogenesis is accelerated dramatically by transgenic co-expression of TAL1 and the catalytic subunit of casein kinase IIHLH transcription factors.

Application Notes

Optimal dilution of the TAL1 antibody should be determined by the researcher.

Immunogen

A recombinant human full-length TAL1 protein was used as the immunogen for the TAL1 antibody.

Storage

Store the TAL1 antibody at 2-8oC (with azide) or aliquot and store at -20oC or colder (without azide).