


Tankyrase 2 Antibody (R35581)

Catalog No.	Formulation	Size
R35581-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	80351
Applications	Western Blot : 1-3ug/ml
Limitations	This Tankyrase 2 antibody is available for research use only.



Western blot testing of transfected HEK293 lysate with Tankyrase 2 antibody at 0.5ug/ml. Lane 1: mock TANK2; Lane 2: TANK1; Lane 3: empty; Lane 4: TANK2. Lower panel shows the same lysates probed for Actin (loading control).

Description

Additional name(s) for this target protein: TRF1-interacting ankyrin-related ADP-ribose polymerase 2, TNKS2

Application Notes

Optimal dilution of the Tankyrase 2 antibody should be determined by the researcher.

1. This Tankyrase 2 antibody will not detect Tankyrase 1.

Immunogen

A portion of amino acids 800-900 was used as the immunogen for this Tankyrase 2 antibody.

Storage

Aliquot and store the Tankyrase 2 antibody at -20oC.