


TBPL1 Antibody (R35312)

Catalog No.	Formulation	Size
R35312-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9519
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This TBPL1 antibody is available for research use only.



Western blot testing of Jurkat lysate with TBPL1 antibody at 2ug/ml. Predicted molecular weight: ~20kDa.

Description

Additional name(s) for this target protein: TBP-like 1, TRF2, TATA box binding protein-related factor 2, TLF

Application Notes

Optimal dilution of the TBPL1 antibody should be determined by the researcher.

Immunogen

Amino acids QIYPFVFESRKEIL were used as the immunogen for this TBPL1 antibody.

Storage

Aliquot and store the TBPL1 antibody at -20oC.