


TBR2 Antibody / EOMES / Eomesodermin (F47508)

Catalog No.	Formulation	Size
F47508-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F47508-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	O95936
Applications	Western Blot : 1:1000
Limitations	This TBR2 antibody is available for research use only.



Description

T-box Brain protein 2, also called EOMES, is a member of a conserved protein family that shares a common DNA-binding domain, the T-box. T-box genes encode transcription factors involved in the regulation of developmental processes. A similar gene disrupted in mice is shown to be essential during trophoblast development and gastrulation.

Application Notes

Titration of the TBR2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 28-57 from the human protein was used as the immunogen for this EOMES / TBR2 antibody.

Storage

Aliquot the TBR2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.