

TCF4 Antibody (R34798)

Catalog No.	Formulation	Size
R34798-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6925
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This TCF4 antibody is available for research use only.



Western blot testing of human heart lysate with TCF4 antibody at 0.5ug/ml. Predicted molecular weight ~71kDa.

Description

Additional name(s) for this target protein: Transcription factor 4, E2-2, ITF2

Application Notes

Optimal dilution of the TCF4 antibody should be determined by the researcher.

1. This TCF4 antibody will detect isoforms a and b.

Immunogen

Amino acids NGGHPSPSRNYGDGT were used as the immunogen for this TCF4 antibody.

Storage

Aliquot and store the TCF4 antibody at -20oC.