


TCF7 Antibody (F40111)

Catalog No.	Formulation	Size
F40111-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40111-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	P36402
Applications	Western Blot : 1:1000
Limitations	This TCF7 antibody is available for research use only.



Description

The T cell specific transcription factor TCF7 activates genes involved in immune regulation and thymocyte differentiation, and is a candidate locus for genetic susceptibility to type 1 diabetes.

Application Notes

Titration of the TCF7 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 3-36 from the human protein was used as the immunogen for this TCF7 antibody.

Storage

Aliquot the TCF7 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.