


## TdT Antibody / Terminal deoxynucleotidyl transferase [clone RM379] (R20398)

Catalog No.	Formulation	Size
R20398-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM379
Purity	Protein A purified from animal origin-free supernatant
UniProt	P04053
Localization	Nuclear
Applications	Immunohistochemistry (FFPE) : 1:100-1:200 Western Blot : 1:1000-1:2000
Limitations	This recombinant TdT antibody is available for research use only.



Western blot testing of human Jurkat cell lysate with recombinant TdT antibody at 1:2000 dilution. Predicted molecular weight: ~56 kDa.



IHC staining of FFPE human thymus tissue with recombinant TdT antibody at 1:200.

## Description

This antibody reacts to human DNA nucleotidyltransferase / Terminal deoxynucleotidyl transferase / TdT.

## Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant TdT antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A peptide corresponding to human Terminal deoxynucleotidyl transferase (TdT) was used as the immunogen for the recombinant TdT antibody.

## Storage

Store the recombinant TdT antibody at -20°C.