


TERT Antibody (R35362)

Catalog No.	Formulation	Size
R35362-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7015
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TERT antibody is available for research use only.



Western blot testing of human skeletal muscle lysate with TERT antibody at 0.5ug/ml. Predicted molecular weight ~126kDa, visualized here at ~140kDa.

Description

Additional name(s) for this target protein: Telomerase reverse transcriptase

Application Notes

Optimal dilution of the TERT antibody should be determined by the researcher.

Immunogen

Amino acids QLRELSEAEVRQHRE were used as the immunogen for this TERT antibody.

Storage

Aliquot and store the TERT antibody at -20oC.