


TFAP2C Antibody / AP2 gamma (R34852)

Catalog No.	Formulation	Size
R34852-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7022
Applications	Immunohistochemistry (FFPE) : 2-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This TFAP2C antibody is available for research use only.



IHC testing of FFPE human skin with TFAP2C antibody at 4ug/ml. HIER: steamed with pH9 Tris/EDTA buffer, HRP-staining.

Description

Additional name(s) for this target protein: Transcription factor AP-2 gamma (activating enhancer-binding protein 2 gamma), AP-2 gamma, TFAP2G

Application Notes

Optimal dilution of the TFAP2C antibody should be determined by the researcher.

Immunogen

Amino acids NKTLEKMEKHRK were used as the immunogen for this TFAP2C antibody.

Storage

Aliquot and store the TFAP2C antibody at -20oC.