


TFF1 Antibody / Estrogen Inducible Protein pS2 [clone ADIB-20] (RQ5039)

Catalog No.	Formulation	Size
RQ5039	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	ADIB-20
Purity	Affinity purified
UniProt	P04155
Applications	Western Blot : 1:500-1:2000
Limitations	This TFF1 antibody is available for research use only.



Western blot testing of human MCF7 cell lysate with TFF1 antibody. Predicted molecular weight ~12 kDa.

Description

Stabilizer of the mucous gel overlying the gastrointestinal mucosa that provides a physical barrier against various noxious agents. May inhibit the growth of calcium oxalate crystals in urine. [UniProt]

Application Notes

Optimal dilution of the TFF1 antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human Estrogen Inducible Protein pS2 / TFF1 was used as the immunogen for the TFF1

antibody.

Storage

Store the TFF1 antibody at -20oC.