


TGFB2 Antibody (pre-pro region) (R36470)

Catalog No.	Formulation	Size
R36470-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Dog, Cow, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	P61812
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TGFB2 antibody is available for research use only.



Description

Additional name(s) for this target protein: Transforming growth factor beta 2; prepro-transforming growth factor beta-2; G-TSF; Glioblastoma-derived T-cell suppressor factor

Application Notes

Optimal dilution of the TGFB2 antibody should be determined by the researcher.

Immunogen

Amino acids STYTSGDQKTIKSTR were used as the immunogen for this TGFB2 antibody. This sequence is from the

latency associated peptide region.

Storage

Aliquot and store the TGFB2 antibody at -20oC.