


TGFBR1 Antibody (R32084)

Catalog No.	Formulation	Size
R32084	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	P36897
Applications	Western Blot : 0.1-0.5ug/ml
Limitations	This TGFBR1 antibody is available for research use only.



Western blot testing of human 1) HeLa and 2) SW579 cell lysate with TGFBR1 antibody. Expected molecular weight ~55 kDa.

Description

Transforming growth factor, beta receptor I is a TGF beta receptor. TGFBR1 is its human gene. The protein encoded by this gene forms a heteromeric complex with type II TGF-beta receptors when bound to TGF-beta, transducing the TGF-beta signal from the cell surface to the cytoplasm. Mutations in this gene have been associated with Loeys-Dietz aortic aneurysm syndrome (LDAS). TGFBR1 regulates cell cycle progression by a unique signaling mechanism that involves its binding to TGFBR2 and activation of TGFBR1. Both are transmembrane serine/threonine receptor kinases. The TGFBR1 receptor may be inactivated in many of the cases of human tumor cells refractory to TGF-beta-mediated cell cycle arrest. Vellucci and Reiss (1997) reported that the TGFBR1 gene is approximately 31 kb long and contains 9 exons. The organization of the segment of the gene that encodes the C-terminal portion of the serine/threonine kinase domain

appears to be highly conserved among members of the gene family.

Application Notes

Optimal dilution of the TGFBR1 antibody should be determined by the researcher.

Immunogen

Amino acids HNRTVIHHRVPNEEDPSLDRPFISEGTTLKDLIYDMTT of human TGFBR1 were used as the immunogen for the TGFBR1 antibody.

Storage

After reconstitution, the TGFBR1 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.