


TIE1 Antibody / Tyrosine-protein kinase receptor Tie-1 (F43869)


Catalog No.	Formulation	Size
F43869-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F43869-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q06806
Applications	Western Blot : 1:1000
Limitations	This TIE1 antibody is available for research use only.



TIE1 antibody western blot analysis in human K562 lysate.



TIE1 antibody western blot analysis in mouse heart tissue lysate.

Description

TIE1 is a transmembrane tyrosine-protein kinase that may modulate TEK/TIE2 activity and contribute to the regulation of angiogenesis. [UniProt]

Application Notes

Titration of the TIE1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 361-389 from the mouse protein was used as the immunogen for this TIE1 antibody.

Storage

Aliquot the TIE1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.