


TIE1 Antibody / Tyrosine-protein kinase receptor Tie-1 (F50681)

Catalog No.	Formulation	Size
F50681-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F50681-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	P35590
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:10-1:50
Limitations	This TIE1 antibody is available for research use only.



Western blot analysis of TIE1 antibody and mouse bladder tissue lysate.



IHC analysis of FFPE human brain tissue stained with TIE1 antibody.

Description

Protein kinases are enzymes that transfer a phosphate group from a phosphate donor, generally the γ phosphate of ATP, onto an acceptor amino acid in a substrate protein. By this basic mechanism, protein kinases mediate most of the signal transduction in eukaryotic cells, regulating cellular metabolism, transcription, cell cycle progression, cytoskeletal rearrangement and cell movement, apoptosis, and differentiation. With more than 500 gene products, the protein kinase family is one of the largest families of proteins in eukaryotes. The family has been classified in 8 major groups based on sequence comparison of their tyrosine (PTK) or serine/threonine (STK) kinase catalytic domains. The tyrosine kinase (TK) group is mainly involved in the regulation of cell-cell interactions such as differentiation, adhesion, motility and death. There are currently about 90 TK genes sequenced, 58 are of receptor protein TK (e.g. EGFR, EPH, FGFR, PDGFR, TRK, and VEGFR families), and 32 of cytosolic TK (e.g. ABL, FAK, JAK, and SRC families).

Application Notes

Titration of the TIE1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 12-42 from the human protein was used as the immunogen for this TIE1 antibody.

Storage

Aliquot the TIE1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.