


TIP47 Antibody (R36268)

Catalog No.	Formulation	Size
R36268-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10226
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This TIP47 antibody is available for research use only.



Western blot testing of A) HeLa and B) human colon cancer lysate with TIP47 antibody at 0.03ug/ml. Image shows isoforms 1-3 being detected.

Description

Additional name(s) for this target protein: Tail interacting protein 47 kD; perilipin 3; PLIN3

Application Notes

Optimal dilution of the TIP47 antibody should be determined by the researcher.

Immunogen

Amino acids TRGAVQSGVDKTKS were used as the immunogen for this TIP47 antibody. This sequence is common to isoforms 1-3.

Storage

Aliquot and store the TIP47 antibody at -20oC.