


TNFR2 Antibody (R34458)


Catalog No.	Formulation	Size
R34458-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	21938
Localization	Cytoplasmic, membrane
Applications	Western Blot : 0.02-0.1ug/ml Immunohistochemistry (FFPE) : 5ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This TNFR2 antibody is available for research use only.



Western blot testing of rat spinal cord lysate with TNFR2 antibody at 0.2ug/ml. Predicted molecular weight: 50/70~80kDa (unmodified/glycosylated).



IHC testing of FFPE human liver tissue with TNFR2 antibody at 5ug/ml. HIER: steam section in pH6 citrate buffer for 20 min and allow to cool prior to staining.

Description

Additional name(s) for this target protein: Tumor necrosis factor receptor superfamily, member 1B, Tumor necrosis factor receptor 2, p75; Tnfrsf1b

Application Notes

Optimal dilution of the TNFR2 antibody should be determined by the researcher.

Immunogen

Amino acids KHFCNKTSDTV CAD were used as the immunogen for this TNFR2 antibody. The amino acid sequence used as immunogen is from the extracellular portion of the protein.

Storage

Aliquot and store the TNFR2 antibody at -20oC.