


## TNNI3 Antibody (R34692)


Catalog No.	Formulation	Size
R34692-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	P19429
Gene ID	7137
Applications	Western Blot : 0.2-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This TNNI3 antibody is available for research use only.



Western blot testing of human heart lysate with TNNI3 antibody at 0.2ug/ml. Predicted molecular weight: ~24kDa, can also observed at ~28kDa.



Western blot of adult mouse (A), fetal mouse (B), adult rat (C) and adult pig (D) heart lysates with TNNI3 antibody at 0.1ug/ml. Predicted molecular weight: ~24kDa, can also observed at ~28kDa.

## Description

Additional name(s) for this target protein: Troponin I type 3

## Application Notes

Optimal dilution of the TNNI3 antibody should be determined by the researcher.

## Immunogen

Amino acids KVTKNITEIAD were used as the immunogen for this TNNI3 antibody.

## Storage

Aliquot and store the TNNI3 antibody at -20oC.