


TP63 Antibody (R33897)

Catalog No.	Formulation	Size
R33897-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8626
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This TP63 antibody is available for research use only.



Western blot testing of A431 lysate with TP63 antibody at 0.3ug/ml

Description

Additional name(s) for this target protein: Transformation related protein 63; Tumor protein 63; p63

Application Notes

Optimal dilution of the TP63 antibody should be determined by the researcher.

1. This TP63 antibody will detect isoforms d, e, f, g, and h with predicted molecular weights between 45~66kDa.

Immunogen

Amino acids ENNAQTQFSEPQYC were used as the immunogen for this TP63 antibody. This sequence is common to isoforms d, e, f, g, and h.

Storage

Aliquot and store the TP63 antibody at -20oC.