


TPH2 Antibody / Tryptophan Hydroxylase 2 (R34844)

Catalog No.	Formulation	Size
R34844-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Pig, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	121278
Applications	Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This TPH2 antibody is available for research use only.



IHC testing of FFPE human brain stem with TPH2 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Tryptophan hydroxylase 2

Application Notes

Optimal dilution of the TPH2 antibody should be determined by the researcher.

Immunogen

Amino acids NKPNSGKNDDKGNK were used as the immunogen for this TPH2 antibody.

Storage

Aliquot and store the TPH2 antibody at -20oC.