


TPPP Antibody (R33583)


Catalog No.	Formulation	Size
R33583-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	11076
Localization	Cytoplasmic, nuclear
Applications	Western Blot : 0.1-0.3ug/ml Immunohistochemistry (FFPE) : 2-5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TPPP antibody is available for research use only.




Western blot testing of human cerebellum lysate with TPPP antibody at 0.001ug/ml.
Predicted molecular weight: ~24 kDa.



Western blot of mouse (A) and rat (B) brain lysate with TPPP antibody at 0.001ug/ml. Predicted molecular weight: ~24 kDa.



IHC staining of FFPE human cerebral cortex with TPPP antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.



IHC staining of FFPE human cerebellum with (top) and without (bottom) TPPP antibody at 2ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Tubulin polymerization promoting protein

Application Notes

Optimal dilution of the TPPP antibody should be determined by the researcher.

Immunogen

Amino acids DPSKDRAAKRLSLES were used as the immunogen for this TPPP antibody.

Storage

Aliquot and store the TPPP antibody at -20oC.