


Transferrin receptor Antibody (R36267)

Catalog No.	Formulation	Size
R36267-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7037
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Transferrin receptor antibody is available for research use only.



Western blot testing of human breast lysate with Transferrin receptor antibody at 1ug/ml.
Predicted molecular weight 85–95 kDa

Description

Additional name(s) for this target protein: TFRC; TRFR; CD71

Application Notes

Optimal dilution of the Transferrin receptor antibody should be determined by the researcher.

Immunogen

Amino acids TRFSLARQVDGDNH were used as the immunogen for this Transferrin receptor antibody.

Storage

Aliquot and store the Transferrin receptor antibody at -20oC.