


## Transglutaminase 6 Antibody (R33519)

Catalog No.	Formulation	Size
R33519-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	343641
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Transglutaminase 6 antibody is available for research use only.



Western blot testing of rat lung lysate with Transglutaminase 6 antibody at 2ug/ml.  
Predicted molecular weight: ~80kDa.

### Description

Additional name(s) for this target protein: Transglutaminase 6; Tgm6; Tg6

### Application Notes

Optimal dilution of the Transglutaminase 6 antibody should be determined by the researcher.

### Immunogen

Amino acids DKYVDSYGRTLED were used as the immunogen for this Transglutaminase 6 antibody.

## Storage

Aliquot and store the Transglutaminase 6 antibody at -20oC.