


TRAPPC2 Antibody / Trafficking protein particle complex 2 [clone 28T69] (RQ8490)

Catalog No.	Formulation	Size
RQ8490	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ug

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	28T69
Purity	Affinity chromatography
UniProt	P0DI81
Applications	Western Blot : 1:500-1:2000
Limitations	This TRAPPC2 antibody is available for research use only.



Western blot testing of human Raji cell lysate with TRAPPC2 antibody. Predicted molecular weight ~16 kDa.

Description

Trafficking protein particle complex 2 may play a role in vesicular transport from endoplasmic reticulum to Golgi. [UniProt]

Application Notes

Optimal dilution of the TRAPPC2 antibody should be determined by the researcher.

Immunogen

A synthetic peptide from human Trafficking protein particle complex 2 protein was used as the immunogen for the TRAPPC2 antibody.

Storage

Store the TRAPPC2 antibody at -20oC.