


TRIF Antibody (R31002)


Catalog No.	Formulation	Size
R31002	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2.5% BSA and 0.025% sodium azide/thimerosal
UniProt	Q8IUC6
Applications	Western Blot : 0.5-1ug/ml IHC (FFPE) : 0.5-1ug/ml IHC (Frozen) : 0.5-1ug/ml
Limitations	This TRIF antibody is available for research use only.



Western blot with TRIF antibody and Lane 1: Jurkat; 2: HL-60 cell lysate.
Predicted/observed molecular weight: 76/70-150kDa.



IHC-P: TRIF antibody testing of human breast cancer tissue



Western blot with TRIF antibody and recombinant human protein; Lane 1: 5ng; 2: 2.5ng; 3: 1.25ng

Description

TIR Domain-containing adaptor molecule 1, also known as TICAM1 or TRIF, is an adapter in responding to activation of toll-like receptors (TLRs). It mediates the delayed cascade of two TLR-associated signaling cascades, where the other one is dependent upon a MyD88 adapter. By genomic sequence analysis, Oshiumi et al.(2003) mapped the gene to chromosome 19p13.3. By coimmunoprecipitation analysis, Oshiumi et al.(2003) showed that the protein interacts specifically with [TLR3](#), but not with other TLRs. Functional analysis showed that the association of TLR3 and TRIF mediates dsRNA activation of IFNB, through either [NFKB](#), [AP1](#), or [IRF3](#). TRIF activation of NFKB was found to occur predominantly through [IRAK1](#) rather than [IRAK2](#). siRNA blockage of the protein reduced IFNB production in response to dsRNA.

Application Notes

The stated application concentrations are suggested starting points. Titration of the TRIF antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

An amino acid sequence from the C-terminus of human TRIF (NNHMWNQRGSQAPEDKTQEAE) was used as the immunogen for this TRIF antibody.

Storage

After reconstitution, the TRIF antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.