


TRIF Antibody (R34971)

Catalog No.	Formulation	Size
R34971-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	148022
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This TRIF antibody is available for research use only.



Western blot testing of Daudi lysate with TRIF antibody at 2ug/ml. Predicted/observed molecular weight: 76/70-150kDa.

Description

Additional name(s) for this target protein: TIR domain containing adaptor inducing interferon-beta, toll-like receptor adaptor molecule 1; TICAM1

Application Notes

Optimal dilution of the TRIF antibody should be determined by the researcher.

Immunogen

Amino acids PDGATFCEDFQVP were used as the immunogen for this TRIF antibody.

Storage

Aliquot and store the TRIF antibody at -20oC.