

TRIF Antibody / TICAM1 (R32274)

Catalog No.	Formulation	Size
R32274	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2.5% BSA and 0.025% sodium azide
UniProt	Q80UF7
Applications	Western Blot : 0.1-0.5ug/ml
Limitations	This TRIF antibody is available for research use only.



Western blot testing of mouse 1) brain, 2) spleen, 3) testis, and 4) liver lysate with TRIF antibody. Predicted molecular weight: 76 kDa but routinely observed at 70-150 kDa.

Description

TICAM1 (TIR domain containing adaptor molecule 1), also known as TRIF, is an adapter in responding to activation of toll-like receptors (TLRs). It mediates the rather delayed cascade of two TLR-associated signaling cascades, where the other one is dependent upon a MyD88 adapter. By genomic sequence analysis, Oshiumi et al. (2003) mapped the TICAM1 gene to chromosome 19p13.3. By coimmunoprecipitation analysis, Oshiumi et al. (2003) showed that TICAM1 interacts specifically with TLR3, but not with other TLRs. Functional analysis showed that the association of TLR3 and TICAM1 mediates dsRNA activation of IFNB, through NFKB, AP1, or IRF3. TICAM1 activation of NFKB was found to occur predominantly through IRAK1 rather than IRAK2. Small interfering (si)RNA blockage of TICAM1, just upstream of the TIR domain, reduced IFNB production in response to dsRNA.

Application Notes

Optimal dilution of the TRIF antibody should be determined by the researcher.

Immunogen

Amino acids QDTEARVSLESLKMNTVAQLVAHQWADMETTE of mouse TRIF were used as the immunogen for the TRIF antibody.

Storage

After reconstitution, the TRIF antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.