


Trim16 Antibody / Tripartite motif-containing protein 16 (RQ6403)


Catalog No.	Formulation	Size
RQ6403	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request


Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q99PP9
Localization	Cytoplasmic
Applications	Western Blot : 0.5-1ug/ml Immunohistochemistry (FFPE) : 2-5ug/ml Immunofluorescence (FFPE) : 5ug/ml Direct ELISA : 0.1-0.5ug/ml
Limitations	This Trim16 antibody is available for research use only.




Immunofluorescent staining of FFPE mouse HEPA1-6 cells with Trim16 antibody (green) and DAPI nuclear stain (blue). HIER: steam section in pH6 citrate buffer for 20 min.



IHC staining of FFPE mouse intestine tissue with Trim16 antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



IHC staining of FFPE rat intestine tissue with Trim16 antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



Western blot testing of 1) rat testis, 2) rat C6, 3) mouse NIH 3T3 and 4) mouse Neuro-2a cell lysate with Trim16 antibody. Predicted molecular weight: 54-63 kDa (multiple isoforms).

Description

Tripartite motif-containing protein 16 is a protein that in humans is encoded by the TRIM16 gene. This gene was identified as an estrogen and anti-estrogen regulated gene in epithelial cells stably expressing estrogen receptor. The protein encoded by this gene contains two B box domains and a coiled-coiled region that are characteristic of the B box zinc finger protein family. The proteins of this family have been reported to be involved in a variety of biological processes including cell growth, differentiation and pathogenesis. Expression of this gene was detected in most tissues. Its function, however, has not yet been determined.

Application Notes

Optimal dilution of the Trim16 antibody should be determined by the researcher.

Immunogen

Recombinant mouse protein (amino acids D65-K210) was used as the immunogen for the Trim16 antibody.

Storage

After reconstitution, the Trim16 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.

