


TRIM7 Antibody (R35198)


Catalog No.	Formulation	Size
R35198-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	81786
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TRIM7 antibody is available for research use only.



Western blot testing of human skeletal muscle lysate with TRIM7 antibody at 1ug/ml.
Predicted molecular weight: ~25kDa.



IHC staining of FFPE human skeletal muscle with TRIM7 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Tripartite motif-containing 7, GNIP1, glycogenin-interacting protein

Application Notes

Optimal dilution of the TRIM7 antibody should be determined by the researcher.

Immunogen

Amino acids QAPAGPPWDITEA were used as the immunogen for this TRIM7 antibody.

Storage

Aliquot and store the TRIM7 antibody at -20oC.