


TRIM71 Antibody (R35741)


Catalog No.	Formulation	Size
R35741-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat, Dog, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	131405
Localization	Cytoplasmic
Applications	Western Blot : 0.3-1ug/ml IHC (FFPE) : 3.75ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This TRIM71 antibody is available for research use only.



Western blot testing of human HEK293 [A] and mouse NIH3T3 [B] lysate with TRIM71 antibody at 0.3ug/ml. Predicted molecular weight: ~93 kDa.



IHC testing of FFPE human cerebellum tissue with TRIM71 antibody at 3.75ug/ml.
Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Additional name(s) for this target protein: Tripartite motif-containing 71

Application Notes

Optimal dilution of the TRIM71 antibody should be determined by the researcher.

Immunogen

Amino acids RQVKAKSLYLQVE were used as the immunogen for this TRIM71 antibody.

Storage

Aliquot and store the TRIM71 antibody at -20oC.