


## TrkA/B/C Antibody / Pan TRK [clone HHH-14] (RQ5057)

| Catalog No. | Formulation                                                                | Size   |
|-------------|----------------------------------------------------------------------------|--------|
| RQ5057      | Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA | 100 ul |

[Bulk quote request](#)

|                    |                                                            |
|--------------------|------------------------------------------------------------|
| Availability       | 1-2 weeks                                                  |
| Species Reactivity | Human                                                      |
| Format             | Purified                                                   |
| Clonality          | Rabbit Monoclonal                                          |
| Isotype            | Rabbit IgG                                                 |
| Clone Name         | HHH-14                                                     |
| Purity             | Affinity purified                                          |
| UniProt            | P04629, Q16620, Q16288                                     |
| Applications       | Western Blot : 1:1000-1:2000                               |
| Limitations        | This TrkA/B/C antibody is available for research use only. |



Western blot testing of human fetal brain lysate with TrkA/B/C antibody. Expected molecular weight: 85~140 kDa depending on glycosylation level.

### Description

TrkA, TrkB and TrkC are members of the neurotrophic tyrosine kinase receptor (NTRK) family.

### Application Notes

Optimal dilution of the TrkA/B/C antibody should be determined by the researcher.

### Immunogen

A synthetic peptide specific to human TrkA/B/C was used as the immunogen for the TrkA/B/C antibody.

## Storage

Store the TrkA/B/C antibody at -20oC.