

## TRPC4AP Antibody (R34529)

Catalog No.	Formulation	Size
R34529-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	26133
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This TRPC4AP antibody is available for research use only.



Western blot testing of U937 lysate with TRPC4AP antibody at 1ug/ml. Predicted molecular weight: ~90kDa.

## Description

Additional name(s) for this target protein: Transient receptor potential cation channel, subfamily C, member 4

## Application Notes

Optimal dilution of the TRPC4AP antibody should be determined by the researcher.

## Immunogen

Amino acids KLTQETTYPNTYIFD were used as the immunogen for this TRPC4AP antibody.

## Storage

Aliquot and store the TRPC4AP antibody at -20oC.