


TRPM7 Antibody (R36190)

Catalog No.	Formulation	Size
R36190-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Dog, Human, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	54822
Applications	Immunohistochemistry (FFPE) : 4-6ug/ml
Limitations	This TRPM7 antibody is available for research use only.



IHC testing of FFPE mouse hippocampus with TRPM7 antibody at 4ug/ml. HIER: steamed with pH6 or pH9 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Transient receptor potential cation channel, subfamily M, member 7

Application Notes

Optimal dilution of the TRPM7 antibody should be determined by the researcher.

Immunogen

Amino acids TKESESTNSVRLML were used as the immunogen for this TRPM7 antibody.

Storage

Aliquot and store the TRPM7 antibody at -20oC.