

Trpv2 Antibody (R35870)

Catalog No.	Formulation	Size
R35870-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	22368
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Trpv2 antibody is available for research use only.



Western blot testing of rat brain lysate with Trpv2 antibody at 2ug/ml. Expected molecular weight: 85-95 kDa depending on glycosylation level.

Description

Additional name(s) for this target protein: Transient receptor potential cation channel, subfamily V, member 2

Application Notes

Optimal dilution of the Trpv2 antibody should be determined by the researcher.

Immunogen

Amino acids SEEDHLPLQVLQSH were used as the immunogen for this Trpv2 antibody.

Storage

Aliquot and store the Trpv2 antibody at -20oC.

References (1)