


TRX Antibody / TXN / Thioredoxin [clone 3C1-H10] (F54050)

Catalog No.	Formulation	Size
F54050-0.1ML	In PBS with 50% glycerol and 0.03% ProClin 300	0.1 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	NA
Format	Purified
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1
Clone Name	3C1-H10
Purity	Protein G affinity
UniProt	P10599
Applications	Western Blot : 1:1000
Limitations	This TRX antibody is available for research use only.



Western blot testing of bacterial lysates containing TRX-tagged proteins of different MWs using TRX antibody at 1:1000.

Description

Application Notes

The stated application concentrations are suggested starting points. Titration of the TRX antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A recombinant protein was used as the immunogen for this TRX antibody.

Storage

Store the TRX antibody at -20oC.