


TSC22D1 Antibody / TSC22 (F53897)

Catalog No.	Formulation	Size
F53897-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F53897-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q15714
Applications	Western Blot : 1:1000-2000
Limitations	This TSC22D1 antibody is available for research use only.



Western blot testing of Jurkat cell lysate with TSC22D1 antibody at 1:2000. Predicted molecular weight: 110 kDa.

Description

TSC22 domain family protein 1 is a transcriptional repressor, acting on the C-type natriuretic peptide (CNP) promoter.

Application Notes

Titration of the TSC22D1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 888-920 from the human protein was used as the immunogen for the TSC22D1 antibody.

Storage

Aliquot the TSC22D1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.