


TTC8 Antibody / BBS8 (R34427)

Catalog No.	Formulation	Size
R34427-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Dog, Pig, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	123016
Applications	Western Blot : 0.2-0.6ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This TTC8 antibody is available for research use only.



Western blot testing of human testis lysate with TTC8 antibody at 0.2ug/ml. Predicted molecular weight: ~58 kDa.

Description

Additional name(s) for this target protein: Tetratricopeptide repeat domain 8; Bardet-Biedl syndrome 8 protein

Application Notes

Optimal dilution of the TTC8 antibody should be determined by the researcher.

Immunogen

Amino acids KEVLKQDNTHVE were used as the immunogen for this TTC8 antibody.

Storage

Aliquot and store the TTC8 antibody at -20oC.