


TTP Antibody / Tristetraprolin / ZFP36 (F51446)

Catalog No.	Formulation	Size
F51446-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F51446-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P26651
Applications	Western Blot : 1:1000
Limitations	This TTP antibody is available for research use only.



Description

Tristetraprolin (TTP), or ZFP36, is a probable regulatory protein with a novel zinc finger structure involved in regulating the response to growth factors. It has been experimentally shown to be able to bind zinc.

Application Notes

Titration of the TTP antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 116-142 from the human protein was used as the immunogen for this Tristetraprolin / TTP antibody.

Storage

Aliquot the TTP antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.