


## TTR Antibody / Transthyretin / Prealbumin (R33909)

Catalog No.	Formulation	Size
R33909-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7276
Localization	Cytoplasmic, membrane
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This TTR antibody is available for research use only.



Western blot testing of human albumin-depleted serum with TTR antibody at 0.03ug/ml.  
Predicted molecular weight: ~16kDa.

### Description

Additional name(s) for this target protein: Prealbumin, Transthyretin

### Application Notes

Optimal dilution of the TTR antibody should be determined by the researcher.

## Immunogen

Amino acids EPFASGKTSESGE were used as the immunogen for this TTR antibody.

## Storage

Aliquot and store the TTR antibody at -20oC.