


TUBGCP2 Antibody / GCP2 (RQ6813)


Catalog No.	Formulation	Size
RQ6813	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q9BSJ2
Localization	Cytoplasmic
Applications	Western Blot : 1-2ug/ml Immunofluorescence (FFPE) : 5ug/ml Direct ELISA : 0.1-0.5ug/ml
Limitations	This TUBGCP2 antibody is available for research use only.



Immunofluorescent staining of FFPE human A431 cells with TUBGCP2 antibody (green) and DAPI nuclear stain (blue). HIER: steam section in pH6 citrate buffer for 20 min.



Western blot testing of human 1) 293T and 2) U-87 MG cell lysate with TUBGCP2 antibody. Predicted molecular weight ~103 kDa.

Description

Gamma-tubulin complex component 2 is a protein that in humans is encoded by the TUBGCP2 gene. It is part of the gamma tubulin complex, which required formicrotubule nucleation at the centrosome. Gamma-tubulin complex is necessary for microtubule nucleation at the centrosome. Gamma-tubulin complex is composed of gamma-tubulin, TUBGCP2, TUBGCP3, TUBGCP4, TUBGCP5 and TUBGCP6.

Application Notes

Optimal dilution of the TUBGCP2 antibody should be determined by the researcher.

Immunogen

Recombinant human protein (amino acids Q213-R888) was used as the immunogen for the TUBGCP2 antibody.

Storage

After reconstitution, the TUBGCP2 antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.